


Product Datasheet

5-HT6 Antibody (orb1523061)

Description	5-HT6 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	Between 1 and 50
Target	5-HT6
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	P50406
NCBI	NP_000862.1
Dilution Range	WB - Not recommended. A304-401A has not performed satisfactorily when used for WB of 5-HT6 in crude preparations (e.g. whole cell lysate). This antibody can be used for WB of enriched (e.g. immunoprecipitated) sources of 5-HT6 .; IP - 2 - 10 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human 5-HT6 by western blot...