


Product Datasheet

MRPP1 Antibody (orb1523080)

Description	MRPP1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	Between 353 and 403
Target	MRPP1
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q7L0Y3
NCBI	NP_060289.2
Dilution Range	WB - Not recommended. orb1523080 has not performed satisfactorily when used for WB of MRPP1 in crude preparations (e.g. whole cell lysate). This antibody can be used for WB of enriched (e.g. immunoprecipitated) sources of MRPP1.; IP - 2 - 10 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human MRPP1 by western blot...