


## Product Datasheet

### BRPF3 Antibody (orb1523636)

|                            |                                                                        |
|----------------------------|------------------------------------------------------------------------|
| <b>Description</b>         | BRPF3 Antibody                                                         |
| <b>Species/Host</b>        | Rabbit                                                                 |
| <b>Reactivity</b>          | Human                                                                  |
| <b>Conjugation</b>         | Unconjugated                                                           |
| <b>Tested Applications</b> | WB                                                                     |
| <b>Immunogen</b>           | Between 1 and 50                                                       |
| <b>Target</b>              | BRPF3                                                                  |
| <b>Preservatives</b>       | Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide |
| <b>Form/Appearance</b>     | Whole IgG                                                              |
| <b>Concentration</b>       | 1000 µg/ml                                                             |
| <b>Storage</b>             | 2 - 8°C                                                                |
| <b>Note</b>                | For research use only                                                  |
| <b>Clonality</b>           | Polyclonal                                                             |
| <b>Uniprot ID</b>          | <a href="#">Q9ULD4</a>                                                 |
| <b>NCBI</b>                | <a href="#">NP_056510.2</a>                                            |
| <b>Dilution Range</b>      | WB - 1:1,000 - 1:5,000                                                 |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                        |



Detection of human BRPF3 by western blot...