


Product Datasheet

ZZEF1 Antibody (orb1523661)

Description	ZZEF1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	Between 525 and 575
Target	ZZEF1
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	O43149
NCBI	NP_055928.3
Dilution Range	WB - Not recommended. A304-068A has not performed satisfactorily when used for WB of ZZEF1 in crude preparations (e.g. whole cell lysate). This antibody can be used for WB of enriched (e.g. immunoprecipitated) sources of ZZEF1.; IP - 2 - 10 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human ZZEF1 by western blot...