


Product Datasheet

TSC22D2 Antibody (orb1523663)

Description	TSC22D2 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	Between 730 and 780
Target	TSC22D2
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	O75157
NCBI	NP_055594.1
Dilution Range	WB - Not recommended. A304-067A has not performed satisfactorily when used for WB of TSC22D2 in crude preparations (e.g. whole cell lysate). This antibody can be used for WB of enriched (e.g. immunoprecipitated) sources of TSC22D2.; IP - 2 - 10 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human TSC22D2 by western blotting