


Product Datasheet

Cyclin M3 Antibody (orb1523667)

Description	Cyclin M3 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	Between 575 and 625
Target	Cyclin M3
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q8NE01
NCBI	NP_060093.3
Dilution Range	WB - Not recommended. Use rabbit anti-Cyclin M3 antibody A304-066A.; IP - 2 - 10 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human Cyclin M3 by western ...