



## Product Datasheet

### PSMD4 Antibody (orb1524052)


|                            |                                                                                                                                                                   |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | PSMD4 Antibody                                                                                                                                                    |
| <b>Species/Host</b>        | Rabbit                                                                                                                                                            |
| <b>Reactivity</b>          | Human                                                                                                                                                             |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                      |
| <b>Tested Applications</b> | IHC, IP, WB                                                                                                                                                       |
| <b>Immunogen</b>           | Between 327 and 377                                                                                                                                               |
| <b>Target</b>              | PSMD4                                                                                                                                                             |
| <b>Preservatives</b>       | Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide                                                                                            |
| <b>Form/Appearance</b>     | Whole IgG                                                                                                                                                         |
| <b>Concentration</b>       | 1000 µg/ml                                                                                                                                                        |
| <b>Storage</b>             | 2 - 8°C                                                                                                                                                           |
| <b>Note</b>                | For research use only                                                                                                                                             |
| <b>Application notes</b>   | Format: Whole IgG                                                                                                                                                 |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                        |
| <b>Uniprot ID</b>          | <b>P55036</b>                                                                                                                                                     |
| <b>NCBI</b>                | <b>NP_002801.1</b>                                                                                                                                                |
| <b>Dilution Range</b>      | WB - 1:2,000 - 1:10,000; IP - 2 - 10 µg/mg lysate; IHC - 1:1,000 - 1:5,000. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE tissue sections. |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                   |



Detection of human PSMD4 by western blot...



Detection of human PSMD4 by western blot...



Detection of human PSMD4 by immunohistoc...