



Product Datasheet

SOD1 Antibody (orb1524136)

Description	SOD1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Between 54 and 104
Target	SOD1
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	P00441
NCBI	NP_000445.1
Dilution Range	WB - 1:2,000 - 1:10,000; IP - Not recommended; IHC - 1:2,000 - 1:10,000. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE tissue sections.
Expiration Date	12 months from date of receipt.



Detection of human SOD1 by western blot....



Detection of human SOD1 by immunohistoch...