



## Product Datasheet

### PDCD2 Antibody (orb1524530)

|                            |                                                                      |
|----------------------------|----------------------------------------------------------------------|
| <b>Description</b>         | PDCD2 Antibody                                                       |
| <b>Species/Host</b>        | Rabbit                                                               |
| <b>Reactivity</b>          | Human, Mouse                                                         |
| <b>Conjugation</b>         | Unconjugated                                                         |
| <b>Tested Applications</b> | IP, WB                                                               |
| <b>Immunogen</b>           | Between 294 and 344                                                  |
| <b>Target</b>              | PDCD2                                                                |
| <b>Preservatives</b>       | Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide |
| <b>Form/Appearance</b>     | Whole IgG                                                            |
| <b>Concentration</b>       | 200 µg/ml                                                            |
| <b>Storage</b>             | 2 - 8°C                                                              |
| <b>Note</b>                | For research use only                                                |
| <b>Application notes</b>   | Format: Whole IgG                                                    |
| <b>Clonality</b>           | Polyclonal                                                           |
| <b>Uniprot ID</b>          | <a href="#">Q16342</a>                                               |
| <b>NCBI</b>                | <a href="#">NP_002589.2</a>                                          |
| <b>Dilution Range</b>      | WB - 1:2,000 - 1:10,000; IP - 2 - 5 µg/mg lysate                     |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                      |



Detection of human and mouse PDCD2 by western blot...



Detection of human PDCD2 by western blot...