


Product Datasheet

IRF3 Antibody (orb1524929)

Description	IRF3 Antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	between 375 and 419
Target	IRF3
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	P70671
NCBI	NP_058545.1
Dilution Range	WB - 1:1,000 - 1:2,500; IP - 2 - 10 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of mouse IRF3 by western blot ...