


Product Datasheet

RREB1 Antibody (orb1525410)

Description	RREB1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	between 1375 and 1425
Target	RREB1
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q92766
NCBI	NP_002946.4
Dilution Range	WB - Not recommended; IP - 2 - 10 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human RREB1 by western blot...