



Product Datasheet

NF90/NF110 Antibody (orb1525424)

Description	NF90/NF110 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	between 325 and 375
Target	NF90/NF110
Preservatives	Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	200 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q12906
NCBI	NP_036350.2
Dilution Range	WB - 1:2,000 - 1:10,000; IP - Not recommended; IHC - 1:100 - 1:1,000. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE tissue sections.
Expiration Date	12 months from date of receipt.



Detection of human and mouse NF90 and NF110 by immuno-blotting.



Detection of mouse NF90 and NF110 by immunohistochemistry.