


Product Datasheet

p150glued Antibody (orb1525510)

Description	p150glued Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	between 1100 and 1150
Target	p150glued
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q14203
NCBI	NP_004073.2
Dilution Range	WB - Not recommended; IP - 2 - 5 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human p150glued by western ...