



Product Datasheet

RNF138 Antibody (orb1525635)


Description	RNF138 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, IP
Immunogen	between 1 and 50
Target	RNF138
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q8WVD3
NCBI	NP_057355.2
Dilution Range	WB - Not recommended. A302-993A has not consistently performed satisfactorily when used for WB of RNF138 in crude preparations (e.g. whole cell lysate). This antibody can be used for WB of enriched (e.g. immunoprecipitated) sources of RNF138.; IP - 2 - 5 µg/mg lysate; IHC - 1:1,000 - 1:5,000. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE tissue sections.
Expiration Date	12 months from date of receipt.



Detection of human RNF138 by western blo...



Detection of mouse RNF138 by immunohisto...



Detection of human RNF138 by immunohisto...