



Product Datasheet

eIF4AIII/EIF4A3 Antibody (orb1525659)


Description	eIF4AIII/EIF4A3 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	between 125 and 175
Target	eIF4AIII/EIF4A3
Preservatives	Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	200 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	P38919
NCBI	NP_055555.1
Dilution Range	WB - 1:2,000 - 1:10,000; IP - 2 - 5 µg/mg lysate; IHC - 1:100 - 1:500. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE tissue sections.
Expiration Date	12 months from date of receipt.



Detection of human and mouse eIF4AIII/EIF4A3



Detection of human eIF4AIII/EIF4A3 by Western blot



Detection of human eIF4AIII/EIF4A3 by immunohistochemistry