



Product Datasheet

eIF4AIII/EIF4A3 Antibody (orb1525660)


Description	eIF4AIII/EIF4A3 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, IP
Immunogen	between 1 and 50
Target	eIF4AIII/EIF4A3
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	P38919
NCBI	NP_055555.1
Dilution Range	WB - Not recommended; IP - 2 - 5 µg/mg lysate; IHC - 1:1,000 - 1:5,000. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE tissue sections.
Expiration Date	12 months from date of receipt.



Detection of human eIF4AIII/EIF4A3 by we...



Detection of mouse eIF4AIII/EIF4A3 by im...



Detection of human eIF4AIII/EIF4A3 by im...