


## Product Datasheet

### L3MBTL3 Antibody (orb1525879)

|                            |                                                                        |
|----------------------------|------------------------------------------------------------------------|
| <b>Description</b>         | L3MBTL3 Antibody                                                       |
| <b>Species/Host</b>        | Rabbit                                                                 |
| <b>Reactivity</b>          | Human                                                                  |
| <b>Conjugation</b>         | Unconjugated                                                           |
| <b>Tested Applications</b> | IP, WB                                                                 |
| <b>Immunogen</b>           | between 1 and 50                                                       |
| <b>Target</b>              | L3MBTL3                                                                |
| <b>Preservatives</b>       | Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide |
| <b>Form/Appearance</b>     | Whole IgG                                                              |
| <b>Concentration</b>       | 1000 µg/ml                                                             |
| <b>Storage</b>             | 2 - 8°C                                                                |
| <b>Note</b>                | For research use only                                                  |
| <b>Application notes</b>   | Format: Whole IgG                                                      |
| <b>Clonality</b>           | Polyclonal                                                             |
| <b>Uniprot ID</b>          | <a href="#">Q96JM7</a>                                                 |
| <b>NCBI</b>                | <a href="#">NP_115814.1</a>                                            |
| <b>Dilution Range</b>      | WB - 1:2,000 - 1:10,000; IP - 2 - 5 µg/mg lysate                       |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                        |



Detection of human L3MBTL3 by western blotting.