


Product Datasheet

MgcRacGAP Antibody (orb1525979)

Description	MgcRacGAP Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	between 150 and 200
Target	MgcRacGAP
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q9H0H5
NCBI	NP_001119575.1
Dilution Range	WB - 1:2,000 - 1:10,000; IP - 5 - 15 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human MgcRacGAP by western ...