


Product Datasheet

BEGAIN Antibody (orb1525989)

Description	BEGAIN Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	between 493 and 543
Target	BEGAIN
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q9BUH8
NCBI	NP_065887.1
Dilution Range	WB - Not recommended; IP - 5 - 15 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human BEGAIN by western blotting