



Product Datasheet

Pygo2 Antibody (orb1526097)

Description	Pygo2 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IP
Immunogen	between 356 and 406
Target	Pygo2
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q9BRQ0
NCBI	NP_612157.1
Dilution Range	WB - Not recommended. Use rabbit anti-Pygo2 antibody A302-730A.; IP - 2 - 5 µg/mg lysate; IHC - 1:500 - 1:2,000. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE tissue sections.
Expiration Date	12 months from date of receipt.



Detection of human Pygo2 by western blot...



Detection of human Pygo2 by immunohistoc...