



## Product Datasheet

### Pygo2 Antibody (orb1526098)

|                            |                                                                                                                                                                                                        |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Pygo2 Antibody                                                                                                                                                                                         |
| <b>Species/Host</b>        | Rabbit                                                                                                                                                                                                 |
| <b>Reactivity</b>          | Human                                                                                                                                                                                                  |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                                                           |
| <b>Tested Applications</b> | IHC, IP                                                                                                                                                                                                |
| <b>Immunogen</b>           | between 356 and 406                                                                                                                                                                                    |
| <b>Target</b>              | Pygo2                                                                                                                                                                                                  |
| <b>Preservatives</b>       | Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide                                                                                                                                 |
| <b>Form/Appearance</b>     | Whole IgG                                                                                                                                                                                              |
| <b>Concentration</b>       | 1000 µg/ml                                                                                                                                                                                             |
| <b>Storage</b>             | 2 - 8°C                                                                                                                                                                                                |
| <b>Note</b>                | For research use only                                                                                                                                                                                  |
| <b>Application notes</b>   | Format: Whole IgG                                                                                                                                                                                      |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                                                             |
| <b>Uniprot ID</b>          | <b>Q9BRQ0</b>                                                                                                                                                                                          |
| <b>NCBI</b>                | <b>NP_612157.1</b>                                                                                                                                                                                     |
| <b>Dilution Range</b>      | WB - Not recommended. Use rabbit anti-Pygo2 antibody A302-730A.; IP - 2 - 5 µg/mg lysate; IHC - 1:500 - 1:2,000. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE tissue sections. |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                                                        |



Detection of human Pygo2 by western blot...



Detection of human Pygo2 by immunohistoc...