



Product Datasheet

UCKL1 Antibody (orb1526220)

Description	UCKL1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IP
Immunogen	between 50 and 100
Target	UCKL1
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q9NWZ5
NCBI	NP_060329.2
Dilution Range	WB - Not recommended. Use rabbit anti-UCKL1 antibody A302-664A.; IP - 2 - 5 µg/mg lysate; IHC - 1:1,000 - 1:5,000. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE tissue sections.
Expiration Date	12 months from date of receipt.



Detection of human UCKL1 by western blot...



Detection of human UCKL1 by immunohistoc...