



Product Datasheet

Pygo1 Antibody (orb1526251)

Description	Pygo1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IP
Immunogen	between 300 and 350
Target	Pygo1
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q9Y3Y4
NCBI	NP_056432.1
Dilution Range	WB - Not recommended. Use rabbit anti-Pygo1 antibody A302-649A.; IP - 2 - 5 µg/mg lysate; IHC - 1:500 - 1:2,000. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE tissue sections.
Expiration Date	12 months from date of receipt.



Detection of human Pygo1 by western blot...



Detection of human Pygo1 by immunohistoc...