



Product Datasheet

CNOT2 Antibody (orb1526410)

Description	CNOT2 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IP
Immunogen	between 490 and 540
Target	CNOT2
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q9NZN8
NCBI	NP_055330.1
Dilution Range	WB - Not recommended. orb1526410 has not performed satisfactorily when used for WB of CNOT2 in crude preparations (e.g. whole cell lysate). This antibody can be used for WB of enriched (e.g. immunoprecipitated) sources of CNOT2.; IP - 5 - 15 µg/mg lysate; IHC - 1:500 - 1:2,000. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE tissue sections.
Expiration Date	12 months from date of receipt.



Detection of human CNOT2 by western blot...



Detection of human CNOT2 by immunohistoc...