



Product Datasheet

CBX7 Antibody (orb1526467)

Description	CBX7 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	between 201 and 251
Target	CBX7
Preservatives	Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	200 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	O95931
NCBI	NP_783640.1
Dilution Range	WB - 1:2,000 - 1:10,000; IP - 5 - 15 µg/mg lysate; IHC - 1:100 - 1:500. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE tissue sections.
Expiration Date	12 months from date of receipt.



Detection of human CBX7 by western blot ...



Detection of human CBX7 by immunohistoch...