



Product Datasheet

CBX2 Antibody (orb1526470)

Description	CBX2 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ChIP, IP
Immunogen	between 482 and 532
Target	CBX2
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q14781
NCBI	NP_005180.1
Dilution Range	WB - Not recommended. orb1526470 has not performed satisfactorily when used for WB of CBX2 in crude preparations (e.g. whole cell lysate). This antibody can be used for WB of enriched (e.g. immunoprecipitated) sources of CBX2.; IP - 2 - 10 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human CBX2 by western blot ...



Localization of CBX2 Binding Sites by Ch...