


Product Datasheet

SAP130 Antibody (orb1526527)

Description	SAP130 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	between 998 and 1048
Target	SAP130
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q9H0E3
NCBI	NP_078821.2
Dilution Range	WB - Not recommended. Use rabbit anti-SAP130 antibody A302-490A.; IP - 2 - 5 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human SAP130 by western blotting