

Product Datasheet

NHE3 antibody (orb152656)

Description	Rabbit polyclonal to NHE3 (ATTO-390). Sodium-hydrogen exchanger 3 (NHE3; Slc9a3) is an epithelial tr...
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	ATTO-390
Tested Applications	ICC, IHC, WB
Immunogen	Synthetic peptide mapping to AA 809 to 831 of rat sequence
Target	NHE3
Preservatives	PBS, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	1 µg/ml of SPC-400 was sufficient for detection of HNE3 in 10 µg of rat kidney tissue lysate by colorimetric immunoblot analysis using Goat anti-rabbit IgG:HRP as the secondary antibody.
Clonality	Polyclonal
MW	84kDa
Uniprot ID	P26433
NCBI	NP_036786
Entrez	24784
Dilution Range	WB (1:1000); IHC (1:25)
Expiration Date	12 months from date of receipt.