


Product Datasheet

GMP synthase Antibody (orb1526665)

Description	GMP synthase Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	between 600 and 650
Target	GMP synthase
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	P49915
NCBI	NP_003866.1
Dilution Range	WB - Not recommended. Use rabbit anti-GMP synthase antibody A302-417A.; IP - 2 - 5 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human GMP synthase by western blot