


Product Datasheet

NPM1 Antibody (orb1526691)

Description	NPM1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	between 100 and 150
Target	NPM1
Preservatives	Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	200 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	P06748
NCBI	NP_002511.1
Dilution Range	WB - 1:2,000 - 1:10,000; IP - Not recommended. Use rabbit anti-NPM1 antibody A302-403A or A302-404A.
Expiration Date	12 months from date of receipt.



Detection of human and mouse NPM1 by wes...