

Product Datasheet

NKCC2 antibody (orb152673)

Description

Rabbit polyclonal to NKCC2 (ATTO-390). NKCC2,

co-tra...

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

ATTO-390

Tested

IHC, WB

Applications

Immunogen

Produced against a synthetic peptide mapping to a segment (amino acids 33-55) of the N-terminal tail of rat NKCC2 (antibody designation L224).

Target

NKCC2

Preservatives

PBS, 50% glycerol, 0.09% sodium azide
*Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

1 µg/ml of SPC-401 was sufficient for detection of NKCC2 in 10 µg of rat kidney tissue lysate by colorimetric immunoblot analysis using Goat anti-rabbit IgG:HRP as the secondary antibody.

Clonality

Polyclonal

MW

160kDa

Uniprot ID

[Q5PQV1](#)

NCBI

[NP_001257546.1](#)

Entrez

[25065](#)

Dilution Range

WB (1:1000); IHC (1:25)

Expiration Date

12 months from date of receipt.