


## Product Datasheet

### BRD3 Antibody (orb1526751)

|                            |                                                                                          |
|----------------------------|------------------------------------------------------------------------------------------|
| <b>Description</b>         | BRD3 Antibody                                                                            |
| <b>Species/Host</b>        | Rabbit                                                                                   |
| <b>Reactivity</b>          | Human                                                                                    |
| <b>Conjugation</b>         | Unconjugated                                                                             |
| <b>Tested Applications</b> | IP                                                                                       |
| <b>Immunogen</b>           | between 275 and 325                                                                      |
| <b>Target</b>              | BRD3                                                                                     |
| <b>Preservatives</b>       | Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide                   |
| <b>Form/Appearance</b>     | Whole IgG                                                                                |
| <b>Concentration</b>       | 1000 µg/ml                                                                               |
| <b>Storage</b>             | 2 - 8°C                                                                                  |
| <b>Note</b>                | For research use only                                                                    |
| <b>Application notes</b>   | Format: Whole IgG                                                                        |
| <b>Clonality</b>           | Polyclonal                                                                               |
| <b>Uniprot ID</b>          | <a href="#">Q15059</a>                                                                   |
| <b>NCBI</b>                | <a href="#">NP_031397.1</a>                                                              |
| <b>Dilution Range</b>      | WB - Not recommended. Use rabbit anti-BRD3 antibody A302-368A.; IP - 5 - 10 µg/mg lysate |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                          |



Detection of human BRD3 by western blot ...