


Product Datasheet

p57kip2 Antibody (orb1526785)

Description	p57kip2 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	between 266 and 316
Target	p57kip2
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	P49918
NCBI	NP_000067.1
Dilution Range	WB - Not recommended. A302-347A has not performed satisfactorily when used for WB of p57KIP2 in crude preparations (e.g. whole cell lysate). This antibody can be used for WB of enriched (e.g. immunoprecipitated) sources of p57KIP2.; IP - 2 - 5 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human p57kip2 by western bl...