



Product Datasheet

POLR3D Antibody (orb1526884)

Description	POLR3D Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	between 175 and 225
Target	POLR3D
Preservatives	Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	200 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	P05423
NCBI	NP_001713.2
Dilution Range	WB - 1:2,000 - 1:10,000; IP - 10 µg/mg lysate; IHC - 1:200 - 1:1,000. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE tissue sections.
Expiration Date	12 months from date of receipt.



Detection of human POLR3D by western blotting



Detection of human POLR3D by immunohistochemistry