


Product Datasheet

MLLT1 Antibody (orb1526929)

Description	MLLT1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	between 150 and 200
Target	MLLT1
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q03111
NCBI	NP_005925.2
Dilution Range	WB - Not recommended. Use rabbit anti-MLLT1 antibody A302-268A.; IP - 3 - 5 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human MLLT1 by western blot...