

Product Datasheet

NCC antibody [Out of stock] (orb152693)

Description

Rabbit polyclonal to NCC (ATTO-594). NCC,

found on the ...

Species/Host

Rabbit

Reactivity

Canine, Human, Mouse, Rat

Conjugation

ATTO-594

Tested

ICC, IEM, IHC, WB

Applications
Immunogen

Produced against a synthetic peptide mapping to a segment of rat NCC (amino acids 74-95), N-terminal

Target

NCC

Preservatives

PBS, 50% glycerol, 0.09% sodium azide
*Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

1 µg/ml of SPC-402 was sufficient for detection of NCC3 in 10 µg of rat kidney tissue lysate by colorimetric immunoblot analysis using Goat anti-rabbit IgG:HRP as the secondary antibody.

Clonality

Polyclonal

MW

160kDa

Uniprot ID
[P55018](#), [P55017](#)
NCBI
[NP_062218](#), [NP_000330](#)
Entrez
[54300](#), [6559](#)
Dilution Range

WB (1:1000), IHC (1:200)

Expiration Date

12 months from date of receipt.