



Product Datasheet

SASH1 Antibody (orb1526932)

Description	SASH1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IP
Immunogen	between 750 and 800
Target	SASH1
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	O94885
NCBI	NP_056093.3
Dilution Range	WB - Not recommended. Use rabbit anti-SASH1 antibody A302-265A.; IP - 3 - 5 µg/mg lysate; IHC - 1:500 - 1:2,000. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE tissue sections.; IHC-IF - 1:250 - 1:1,000
Expiration Date	12 months from date of receipt.



Detection of human SASH1 by western blot...



Detection of human SASH1 by immunohistoc...