



Product Datasheet

DOCK4 Antibody (orb1526938)

Description	DOCK4 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	between 1725 and 1775
Target	DOCK4
Preservatives	Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	200 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q8N1I0
NCBI	NP_055520.3
Dilution Range	WB - 1:2,000 - 1:10,000; IP - 3 - 5 µg/mg lysate; IHC - 1:100 - 1:500. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE tissue sections.
Expiration Date	12 months from date of receipt.



Detection of human DOCK4 by western blot...



Detection of human DOCK4 by immunohistoc...