


Product Datasheet

CAMSAP1 Antibody (orb1526945)

Description	CAMSAP1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	between 1050 and 1100
Target	CAMSAP1
Preservatives	Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	200 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q5T5Y3
NCBI	NP_056262.3
Dilution Range	WB - 1:2,000 - 1:10,000; IP - 3 - 5 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human CAMSAP1 by western bl...