


Product Datasheet

ARID4B Antibody (orb1526988)

Description	ARID4B Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	between 450 and 500
Target	ARID4B
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q4LE39
NCBI	NP_057458.4
Dilution Range	WB - Not recommended. Use rabbit anti-ARID4B antibody A302-233A.; IP - 5 - 10 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human ARID4B by western blotting