



Product Datasheet

NOPP140 Antibody (orb1527063)

Description	NOPP140 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IP
Immunogen	between 1 and 50
Target	NOPP140
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q14978
NCBI	NP_004732.2
Dilution Range	WB - Not recommended. Use rabbit anti-NOPP140 antibody A302-185A.; IP - 5 - 15 µg/mg lysate; IHC - 1:500 - 1:2,000. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE tissue sections.
Expiration Date	12 months from date of receipt.



Detection of human NOPP140 by western bl...



Detection of human NOPP140 by immunohist...