



Product Datasheet

IP3R1 Antibody (orb1527112)


Description	IP3R1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	between 2708 and 2758
Target	IP3R1
Preservatives	Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	200 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q14643
NCBI	Q14643.2
Dilution Range	WB - 1:2,000 - 1:10,000; IP - 2 - 5 µg/mg lysate; IHC - 1:100 - 1:500. Epitope retrieval with Tris-EDTA pH 9.0 is recommended for FFPE tissue sections.; IHC-IF - 1:50 - 1:250
Expiration Date	12 months from date of receipt.



Detection of human IP3R1 by western blot...



Detection of human IP3R1 by western blot...



Detection of mouse IP3R1 by immunohistoc...