


Product Datasheet

CKAP2L Antibody (orb1527122)

Description	CKAP2L Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	between 200 and 250
Target	CKAP2L
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q8IYA6
NCBI	NP_689728.3
Dilution Range	WB - Not recommended. orb1527122 has not performed satisfactorily when used for WB of CKAP2L in crude preparations (e.g. whole cell lysate). This antibody can be used for WB of enriched (e.g. immunoprecipitated) sources of CKAP2L.; IP - 2 - 5 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human CKAP2L by western blotting