


Product Datasheet

CKAP2L Antibody (orb1527123)

Description	CKAP2L Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP
Immunogen	between 25 and 75
Target	CKAP2L
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q8IYA6
NCBI	NP_689728.3
Dilution Range	WB - Not recommended; IP - 2 - 5 µg/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human CKAP2L by western blo...