

Product Datasheet

ENaC alpha antibody (orb152717)

Description

Rabbit polyclonal antibody against ENaC

Species/Host	Rabbit
Reactivity	Frog, Human, Mouse, Rat
Conjugation	Biotin
Tested Applications	ICC, IHC, IP, WB
Immunogen	Produced against a synthetic peptide mapping to the N-terminal of the alpha subunit (amino acids 46-68) of rat Alpha ENaC (antibody designation 3560-2).
Target	ENaC
Preservatives	PBS, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	1 µg/ml of SPC-403 was sufficient for detection of alpha-ENaC in 35 µg of rat kidney tissue lysate by colorimetric immunoblot analysis using Goat anti-rabbit IgG:HRP as the secondary antibody.
Clonality	Polyclonal
MW	85kDa
Uniprot ID	Q6IRJ1
NCBI	NP_113736
Entrez	25122
Dilution Range	WB (1:1000), ICC/IF (1:400); IHC (1:25)
Expiration Date	12 months from date of receipt.