



Product Datasheet

QSER1 Antibody (orb1527274)

Description	QSER1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IP
Immunogen	between 1407 and 1457
Target	QSER1
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q2KHR3
NCBI	NP_001070254.1
Dilution Range	WB - Not recommended. Use rabbit anti-QSER1 antibody A302-071A.; IP - 2 - 5 µg/mg lysate; IHC - 1:500 - 1:2,000. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE tissue sections.
Expiration Date	12 months from date of receipt.



Detection of human QSER1 by western blot...



Detection of human QSER1 by immunohistoc...