



Product Datasheet

RALY Antibody (orb1527277)


Description	RALY Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	between 256 and 306
Target	RALY
Preservatives	Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	200 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q9UKM9
NCBI	NP_057951.1
Dilution Range	WB - 1:2,000 - 1:10,000; IP - Not recommended. Use rabbit anti-RALY antibody A302-069A.; IHC - 1:200 - 1:1,000. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE tissue sections.
Expiration Date	12 months from date of receipt.



Detection of human RALY by western blot....



Detection of human RALY by western blot ...



Detection of mouse RALY by immunohistoch...