



## Product Datasheet

### RALY Antibody (orb1527278)


|                            |                                                                                                                                                                                                       |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | RALY Antibody                                                                                                                                                                                         |
| <b>Species/Host</b>        | Rabbit                                                                                                                                                                                                |
| <b>Reactivity</b>          | Human, Mouse                                                                                                                                                                                          |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                                                          |
| <b>Tested Applications</b> | IHC, IP, WB                                                                                                                                                                                           |
| <b>Immunogen</b>           | between 256 and 306                                                                                                                                                                                   |
| <b>Target</b>              | RALY                                                                                                                                                                                                  |
| <b>Preservatives</b>       | Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide                                                                                                                                  |
| <b>Form/Appearance</b>     | Whole IgG                                                                                                                                                                                             |
| <b>Concentration</b>       | 200 µg/ml                                                                                                                                                                                             |
| <b>Storage</b>             | 2 - 8°C                                                                                                                                                                                               |
| <b>Note</b>                | For research use only                                                                                                                                                                                 |
| <b>Application notes</b>   | Format: Whole IgG                                                                                                                                                                                     |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                                                            |
| <b>Uniprot ID</b>          | <b>Q9UKM9</b>                                                                                                                                                                                         |
| <b>NCBI</b>                | <b>NP_057951.1</b>                                                                                                                                                                                    |
| <b>Dilution Range</b>      | WB - 1:2,000 - 1:10,000; IP - Not recommended. Use rabbit anti-RALY antibody A302-069A.; IHC - 1:200 - 1:1,000. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE tissue sections. |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                                                       |



Detection of human RALY by western blot....



Detection of human RALY by western blot ...



Detection of human RALY by immunohistochem...