



## Product Datasheet

### RelA, Phospho (S536) Antibody (orb1527302)

|                            |                                                                                                                                                        |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | RelA, Phospho (S536) Antibody                                                                                                                          |
| <b>Species/Host</b>        | Rabbit                                                                                                                                                 |
| <b>Reactivity</b>          | Human                                                                                                                                                  |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                           |
| <b>Tested Applications</b> | IHC, WB                                                                                                                                                |
| <b>Immunogen</b>           | surrounding serine 536                                                                                                                                 |
| <b>Target</b>              | RelA, Phospho (S536)                                                                                                                                   |
| <b>Preservatives</b>       | Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide                                                                                   |
| <b>Form/Appearance</b>     | Whole IgG                                                                                                                                              |
| <b>Concentration</b>       | 200 µg/ml                                                                                                                                              |
| <b>Storage</b>             | 2 - 8°C                                                                                                                                                |
| <b>Note</b>                | For research use only                                                                                                                                  |
| <b>Application notes</b>   | Format: Whole IgG                                                                                                                                      |
| <b>Clonality</b>           | Polyclonal                                                                                                                                             |
| <b>Uniprot ID</b>          | <a href="#">Q04206</a>                                                                                                                                 |
| <b>NCBI</b>                | <a href="#">NP_068810.2</a>                                                                                                                            |
| <b>Dilution Range</b>      | WB - 1:2,000 - 1:10,000; IP - Not recommended; IHC - 1:200 - 1:1,000. Epitope retrieval with Tris-EDTA pH 9.0 is recommended for FFPE tissue sections. |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                        |



Detection of human phospho RelA-S536 by ...



Detection of human Phospho RelA (S536) b...