



Product Datasheet

Cyclin L1 Antibody (orb1527307)

Description	Cyclin L1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IP
Immunogen	between 350 and 400
Target	Cyclin L1
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q9UK58
NCBI	NP_064703.1
Dilution Range	WB - Not recommended. Use rabbit anti-Cyclin L1 antibody A302-058A.; IP - 10 µg/mg lysate; IHC - 1:500 - 1:2,000. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE tissue sections.
Expiration Date	12 months from date of receipt.



Detection of human Cyclin L1 by western ...



Detection of human Cyclin L1 by immunohi...